

DataHarmonizer: Manual for Emissions and Remediation of Airborne PCBs in Schools Data (Project 3)

Purpose

The ISRP DataHarmonizer is a standardized browser-based spreadsheet editor and validator that can reinforce using standardized columns (AKA fields, variables) and controlled vocabulary for certain values (e.g., NA for missing values). Using the DataHarmonizer can also help ensure your data meets our ISRP tabular data structure guideline and improve data exchange between ISRP research projects and support cores.

Template

The templates were created using datasets from Jason Hua's PCB in emission and remediation of airborne in Project 3. Each dataset should have its own template.

Procedure

1. Extract the zip file and navigate into the extracted folder. Open main.html. The DataHarmonizer will open in your default browser. Click **File** on the top left toolbar, and then click **Load template from menu**. Now you'll see a pop-up window. Pull the downwards arrow to change the template, and then click **Open**.
2. To import data from local xlsx, xls, tsv and csv files to the DataHarmonizer, click **File** on the top left toolbar, and then click **Open**. You can also copy and paste data from a spreadsheet to the DataHarmonizer.
3. Definitions, examples, and guidance about field headers are available by double clicking on the field headers, or by using the Reference Guide. To access the Reference Guide, click on the **Help** button, and then click **Reference Guide**.

Required fields are organized into subsections. See Appendix A for required field definitions and guidance.

Required Fields	Field Name	Subsection	Template Name
sample id	sample_id	Sampling Metadata	Jason_emissions Jason_concentrations Jason_sampleType

sample type	Sample_type	Sampling Metadata	Jason_sampleType
room id	room_id	Sampling Metadata	Jason_sampleType
materials	materials	Sampling Metadata	Jason_sampleType
school id	school_id	Sampling Metadata	Jason_sampleType
sampling round	sampling_round	Sampling Metadata	Jason_sampleType
remediation	remediation	Sampling Metadata	Jason_sampleType
sampling start date	sampling_start_date	Sampling Metadata	Jason_sampleType
sampling end date	sampling_end_date	Sampling Metadata	Jason_sampleType
remediation method	remediation_method	Sampling Metadata	Jason_sampleType
tPCB	tPCB	PCB Measurements	Jason_emissions Jason_concentrations
PCB1-PCB209	PCB1-PCB209	PCB Measurements	Jason_emissions Jason_concentrations

All PCB analyte fields are required. If a value is unavailable, enter NA. If the analyte is below the limit of quantification, enter ND. Do not leave cells blank.

Some fields contain a picklist for predefined value. For example, PCB fields contain ND (under detection) and NA (missing value).

4. Validate the imported data by clicking on the **Validate** button on the top left toolbar. Missing values and invalid entries in required fields will be highlighted in red. Address any invalid data that was flagged in red in the template.
 - Pale Red = Incorrect data format
 - Dark Red = Required data missing

Note: It is possible to export incomplete or invalid data. Make sure to review any errors prior to exporting.
5. Once you've completed validating your data, save data by clicking **File**, and then **Save As**. Enter the file name, pull the downwards arrow to select the file format (.xlsx or .csv), and press **Save**. The file will be automatically saved in Downloads. Files saved in the DataHarmonizer can be re-opened directly.